

Anti-Tumor Necrosis Factor Receptor antibody [H398] (STJ16100922)

STJ16100922

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tumor Necrosis Factor Receptor is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F/FC/ELISA/IP/WB
Host/Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H398
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence