

Anti-TNF-RI antibody [HM104] (STJ16100930)

STJ16100930

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-TNF-RI is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA and Immunoprecipitation research applications. |
| Applications | IHC-F/FC/ELISA/IP |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Monoclonal |
| Clone ID | HM104 |
| Concentration | 100 ug/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.1% BSA and 0.02% sodium azide |
| Isotype | IgG2a |
| Storage | Store at 2-8°C upon receipt. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|------------------------------|
| Gene ID | 7132 |
| Gene Symbol | TNFRSF1A |
| Uniprot ID | TNFR1A_HUMAN |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |