

**Anti-MRP4 antibody [M4I-80] (STJ16100936)**

STJ16100936

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-MRP4 is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
<b>Applications</b>	ICC/IHC-F/WB/FC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human/Mouse/Monkey

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	M4I-80
<b>Concentration</b>	200 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	ICC IHC-F 1:20 WB 1:50 FC
<b>Formulation</b>	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
<b>Isotype</b>	IgG2a
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">10257</a>
<b>Gene Symbol</b>	<a href="#">ABCC4</a>
<b>Uniprot ID</b>	<a href="#">MRP4_HUMAN</a>
<b>Immunogen</b>	MRP4 aa 372-431
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	