

Anti-DR5 antibody [DR5-01] (STJ16100937)

STJ16100937

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DR5 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DR5-01
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	8795
Gene Symbol	TNFRSF10B
Uniprot ID	TR10B_HUMAN
Immunogen	Recombinant fusion protein of human IgG heavy chain and extracellular domain of DR5.
Immunogen Region	
Specificity	
Immunogen Sequence	