

Anti-hRPP20 Immunoglobulin antibody [14B11] (STJ16100943)

STJ16100943

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-hRPP20 Immunoglobulin is suitable for use in ELISA, Western Blot, Immunofluorescence and Immunoprecipitation research applications.
Applications	ELISA/WB/IF/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	14B11
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg BSA
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	GST-hRPP20
Immunogen Region	
Specificity	
Immunogen Sequence	