

**Anti-Collagenase G antibody [1D2] (STJ16100946)**

STJ16100946

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Collagenase G is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	1D2
<b>Concentration</b>	0.75 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	ELISA1:500-1:5000 WB1:2000-1:5000
<b>Formulation</b>	PBS, pH 7.2
<b>Isotype</b>	IgG2b
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Collagenase from Clostridium histolyticum
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	