

**Anti-PTPRR protein isoforms antibody [6A6] (STJ16100949)**

STJ16100949

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-PTPRR protein isoforms is suitable for use in Immunohistochemistry, Western Blot and Immunoprecipitation research applications.
<b>Applications</b>	IHC/WB/IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	6A6
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg BSA
<b>Isotype</b>	IgG2
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">5801</a>
<b>Gene Symbol</b>	<a href="#">PTPRR</a>
<b>Uniprot ID</b>	<a href="#">PTPRR_HUMAN</a>
<b>Immunogen</b>	Recombinant PTP-SL-GST (NCBI accession number BN000437)
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	