

**Anti-Citrullinated Fibrinogen Immunoglobulin antibody [3F2] (STJ16100954)**

STJ16100954

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Citrullinated Fibrinogen Immunoglobulin is suitable for use in Immunohistochemistry, Western Blot, ELISA and Immunoprecipitation research applications.
<b>Applications</b>	IHC/WB/ELISA/IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Mouse/Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	3F2
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	ELISA 50-500 ng/μl
<b>Formulation</b>	PBS
<b>Isotype</b>	IgG2ak
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Deiminated murine fibrinogen peptide
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	