

Anti-Prostate Specific Antigen antibody [ER-PR8] (STJ16100961)

STJ16100961

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Prostate Specific antigen is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-F/IHC-P/IF/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ER-PR8
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P 1:10 IF WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	354
Gene Symbol	KLK3
Uniprot ID	KLK3_HUMAN
Immunogen	Isolated from human prostatic carcinoma
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	