

Anti-MRP8 antibody [M8I-74] (STJ16100963)

STJ16100963

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP8 is suitable for use in Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
Applications	ICC/IHC-F/WB
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M8I-74
Concentration	250 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	ICC IHC-F 1:20 WB 1:50
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6279
Gene Symbol	S100A8
Uniprot ID	S10A8_HUMAN
Immunogen	MRP8 aa 1-83
Immunogen Region	
Specificity	
Immunogen Sequence	