

Anti-MRP7 antibody [M7I-3] (STJ16100967)

STJ16100967

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP7 is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	ICC/IHC-F/WB/FC
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M7I-3
Concentration	250 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	ICC IHC-F 1:20 WB FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	89845
Gene Symbol	ABCC10
Uniprot ID	MRP7_HUMAN
Immunogen	human MRP7 aa 194-272
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	