

Anti-P53 antibody [BP53+12] (STJ16100968)

STJ16100968

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-P53 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-F/IHC-P/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BP53+12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:5-1:10
Formulation	Contains 0.09% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7157
Gene Symbol	TP53
Uniprot ID	P53_HUMAN
Immunogen	Human recombinant p53 protein
Immunogen Region	
Specificity	
Immunogen Sequence	