

## Anti-Lymphocytes antibody [CT-4] (STJ16100971)

## **GENERAL INFORMATION**

Product Type Primary antibodies

Short Description Mouse monoclonal antibody anti-Lymphocytes is suitable for use in Immunohistochemistry and Flow Cytometry research

applications.

Applications IHC-F/FC Host/Source Mouse

Reactivity Guinea Pig/Human/Mouse/Rat

## **PRODUCT PROPERTIES**

Clonality Monoclonal Clone ID CT-4 Concentration 500 ug/mL Conjugation Unconjugated Purification

Dilution Range IHC-F 1:100-1:1000

FC

Formulation cell culture supernatant with 0.7% BSA and 0.1% sodium azide

Isotype IgG3
Storage Instruction Store at 2-8°C upon receipt.

## **TARGET INFORMATION**

Gene ID Gene Symbol Uniprot ID Immunogen Immunogen Region Specificity Immunogen Sequence