

Anti-MRP6 antibody [M6II-68] (STJ16100976)

STJ16100976

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP6 is suitable for use in Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
Applications	ICC/IHC-F/WB
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M6II-68
Concentration	50 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	ICC IHC-F 1:20-1:50 WB
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	368
Gene Symbol	ABCC6
Uniprot ID	MRP6_HUMAN
Immunogen	mouse MRP6 aa 843-999
Immunogen Region	
Specificity	
Immunogen Sequence	