

Anti-Fibronectin antibody [568] (STJ16100986)

STJ16100986

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Fibronectin is suitable for use in Flow Cytometry, Western Blot, Immunofluorescence and Immunohistochemistry research applications.
Applications	FC/WB/IF/IHC
Host/Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	568
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution	Flow cytometry (1-2 µg/million cells in 0, 1 ml). Immunoblotting (1 µg/ml for 2h at RT). Immunofluorescence (1-2 µg/ml).
Range	Immunohistology (1-2 µg/ml for 30-60 min at RT; staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris PBS with 0.02% Sodium Azide.
Formulation	
Isotype	IgG1k
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	2335
Gene Symbol	FN1
Uniprot ID	F1NC_HUMAN
Immunogen	A BALB/c mouse was immunized with a high molecular fraction of proteins secreted by cultured human fibroblasts. Fusion partner: P3-C63-Ag8.653.
Immunogen Region	
Specificity	
Immunogen Sequence	

