

Anti-Progesterone antibody [EBS-O-048] (STJ16100991)

STJ16100991

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Progesterone is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	All Species

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-O-048
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	10424
Gene Symbol	PGRMC2
Uniprot ID	PGRMC2_HUMAN
Immunogen	A BALB/c mouse was immunized with Progesterone-11a-hemisuccinate conjugated to BSA (Progesterone-11a-HMS-BSA). Fusion partner: X63-Ag8-653.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

