

Anti-MHC II DR-alpha + beta chain antibody [Bra30] (STJ16101040)

STJ16101040

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MHC II DR-alpha + beta chain is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Bra30
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunohistology (1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 mi
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2ak
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A BALB/c mouse was immunized with REH (human ALL cell line). Fusion partner: SP-2/0.
Immunogen Region	
Specificity	
Immunogen Sequence	

