

Anti-CD74 antibody [EBS-O-041] (STJ16101044)

STJ16101044

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD74 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	FC/IHC-F/IF/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-O-041
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1K
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	972
Gene Symbol	CD74
Uniprot ID	HG2A_HUMAN
Immunogen	B lymphoblastoid cell line HFB1.
Immunogen Region	
Specificity	
Immunogen Sequence	

