

## Anti-Cytomegalovirus p65 antibody [CMV221] (STJ16101049)

STJ16101049

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Cytomegalovirus p65 is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	ELISA/FC/IHC-F/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	CMV

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	CMV221
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	IgG1K
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Recombinant CMV-p65 protein
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

