

Anti-MFG-E8/Lactadherin/p47 antibody [MFG-06] (STJ16101070)

STJ16101070

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MFG-E8/Lactadherin/p47 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	FC/IHC-F/IHC-P/IF/WB
Host/Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MFG-06
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1K
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4240
Gene Symbol	MFGE8
Uniprot ID	MFGM_HUMAN
Immunogen	milk fat globule membrane preparation derived from human breast cancer cell line PMC42
Immunogen Region	
Specificity	
Immunogen Sequence	

