

Anti-XRCC5/6-Ku-holo; p70/86 antibody [EBS-C-003] (STJ16101072)

STJ16101072

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-XRCC5/6-Ku-holo; p70/86 is suitable for use in Immunofluorescence and Immunohistochemistry research applications.
Applications	IF/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-C-003
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Immunofluorescence (0.5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30-60 minutes at RT).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A BALB/c mouse was immunized with HeLa cells. Fusion partner: SP2/0.
Immunogen Region	
Specificity	
Immunogen Sequence	

