

Anti-Liver Antigen antibody [EP-1] (STJ16101073)

STJ16101073

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Liver antigen is suitable for use in Flow Cytometry, Immunofluorescence and Immunohistochemistry research applications.
Applications	FC/IF/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EP-1
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (0.5-1.0 µg/million cells in 0.1 ml). Immunofluorescence (0.5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; an appropriate antigen retrieval method for staining of formalin-fixed tissues has not been established to date).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2bk
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A BALB/c mouse was immunized with SK-H1A9-2 human hepatocarcinoma cells. Fusion partner: NS-1.
Immunogen Region	
Specificity	
Immunogen Sequence	

