

Anti-NAB2 antibody [EBS-T-005] (STJ16101075)

STJ16101075

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-NGFI-A-binding protein 2 is suitable for use in Immunohistochemistry research applications.
Applications	IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-T-005
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Immunohistology (1-2 µg/ml for 30 minutes at RT).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4665
Gene Symbol	NAB2
Uniprot ID	NAB2_HUMAN
Immunogen	A BALB/c mouse was immunized with partial recombinant NAB2. Fusion partner: Sp2/0.
Immunogen Region	
Specificity	
Immunogen Sequence	

