

Anti-MALT1 antibody [EBS-T-007] (STJ16101077)

STJ16101077

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-MALT1 is suitable for use in Western Blot and Immunohistochemistry research applications. |
| Applications | WB/IHC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | EBS-T-007 |
| Concentration | Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL. |
| Conjugation | Unconjugated |
| Purification | Affinity purified from tissue culture. |
| Dilution Range | Immunoblotting (1-2 µg/ml). Immunohistology (formalin-fixed: 1-2 µg/ml for 30 min at RT; requires boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10-20 min. followed by cooling at RT for 20 min). |
| Formulation | PBS with 0.02% Sodium Azide. |
| Isotype | IgG1k |
| Storage Instruction | Store for up to 1 year at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 10892 |
| Gene Symbol | MALT1 |
| Uniprot ID | MALT1_HUMAN |
| Immunogen | A BALB/c mouse was immunized with human MALT1 recombinant protein. Fusion partner: Sp2/0. |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |

