

Anti-CD56 antibody [EBS-CD-033] (STJ16101115)

STJ16101115

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD56 is suitable for use in ELISA, Flow Cytometry, Western Blot, Immunofluorescence and Immunohistochemistry research applications.
Applications	ELISA/FC/WB/IF/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-CD-033
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution	ELISA (solid phase; not known; tracer: 0, 001-100 µg/ml for 30 min at RT). Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml).
Range	Immunoblotting (1-2 µg/ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (formalin-fixed: 1-2 µg/ml for 30 min
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4684
Gene Symbol	NCAM1
Uniprot ID	NCAM1_HUMAN
Immunogen	A BALB/c mouse was immunized with a crude membrane preparation of a small cell lung cancer. Fusion partner: SP2/0.
Immunogen Region	
Specificity	
Immunogen Sequence	