

**Anti-CD8a antibody [IBL-3/25] (STJ16101130)**

STJ16101130

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-CD8a is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	FC/IHC-F/IF
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	IBL-3/25
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	IgGK
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	925
<b>Gene Symbol</b>	CD8A
<b>Uniprot ID</b>	CD8A_HUMAN
<b>Immunogen</b>	mouse lymph node cells
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	