

Anti-Interleukin-6 antibody (STJ16101142)

STJ16101142

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Interleukin-6 is suitable for use in FUNC, ELISA and Western Blot research applications.
Applications	FUNC/ELISA/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
Uniprot ID	IL6_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	