

**Anti-ICAM-1 antibody {Biotin} (STJ16101145)**

STJ16101145

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-ICAM-1 is suitable for use in Immunohistochemistry, ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	IHC-F/ELISA/IP/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Baboon/Chimpanzee/Cynomolgus Monkey/Rhesus

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	50 ug/mL
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	Ig
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	3383
<b>Gene Symbol</b>	ICAM1
<b>Uniprot ID</b>	ICAM1_HUMAN
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	