

Anti-Intestinal Fatty Acid Binding Protein antibody (STJ16101147)

STJ16101147

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Intestinal Fatty Acid Binding Protein is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat/Sheep/Swine

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081