

Anti-Secretory Leucocyte Protease Inhibitor antibody (STJ16101148)

STJ16101148

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Secretory Leucocyte Protease Inhibitor is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence