

Anti-Desmin K5 antibody (STJ16101153)

STJ16101153

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Desmin K5 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/IHC-P/WB
Host/Source	Rabbit
Reactivity	Human/Bovine/Caprine/Chicken/Hamster/Mouse/Rat/Zebrafish

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	IHC-F 1:25-1:100 IHC-P 1:25-1:100 WB 1:100-1:500
Formulation	PBS with 0.09% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	anti Chicken gizzard muscle desmin
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081