

Anti-BLNK antibody (STJ16101302)

STJ16101302

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BLNK is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	29760
Gene Symbol	BLNK
Uniprot ID	BLNK_HUMAN
Immunogen	A fusion protein representing amino acids 171-356 of human BLNK.
Immunogen Region	
Specificity	
Immunogen Sequence	