

## Anti-C1q antibody {FITC} (STJ16101305)

STJ16101305

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Sheep polyclonal antibody anti-C1q is suitable for use in Immunohistochemistry research applications.
<b>Applications</b>	IHC-F
<b>Host/Source</b>	Sheep
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	10 mg/mL
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	
<b>Formulation</b>	PBS
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Human C1q, purified from plasma.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	