

Anti-Lactoferrin antibody (STJ16101307)

STJ16101307

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse polyclonal antibody anti-Lactoferrin is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4057
Gene Symbol	LTF
Uniprot ID	TRFL_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	