

Anti-caspase-7 antibody (C-Term) (STJ16101313)

STJ16101313

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-caspase-7 (C-Term) is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC-P
Host/Source	Rabbit
Reactivity	

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	840
Gene Symbol	CASP7
Uniprot ID	CASP7_HUMAN
Immunogen	A 16 amino acid peptide from human Caspase-7 carboxy-terminus
Immunogen Region	C-Term
Specificity	
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081