

Anti-PGP 9.5 antibody (STJ16101345)

STJ16101345

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PGP 9.5 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	IHC-P 1:100-1:500
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7345
Gene Symbol	UCHL1
Uniprot ID	UCHL1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

