

Anti-E-Cadherin antibody [BS38] (STJ16101355)

STJ16101355

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-E-Cadherin is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human/Dog/Rat/Mouse/Pig/Sheep/Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BS38
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:100-1:400
Formulation	TRIS with 0.03% sodium azide, pH7, 2
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	999
Gene Symbol	CDH1
Uniprot ID	CADH1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081