

Anti-Pan Cytokeratin antibody [BS5] (STJ16101360)

STJ16101360

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Pan Cytokeratin is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-P/IHC-F
Host/Source	Mouse
Reactivity	Human/Sheep/Pig/Dog/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BS5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:100-1:400 IHC-F
Formulation	TRIS with 0.03% sodium azide, pH7, 2
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence