

Anti-Isocitrate Dehydrogenase 1 IDH1-R132H antibody [HMab-1] (STJ16101369)

STJ16101369

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Isocitrate Dehydrogenase 1 IDH1-R132H is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HMab-1
Concentration	Lot specific
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:25-1:50
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	KLH-conjugated peptide human IDH1
Immunogen Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081