

Anti-HCL antibody [DBA.44] (STJ16101378)

STJ16101378

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Hairy Cell Leukemia is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DBA.44
Concentration	Lot specific
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:15-1:30
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgMk
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	DEAU cell line established from a large cell lymphoma
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081