

Anti-CD33 antibody [WM53] {PE} (STJ16101385)

STJ16101385

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD33 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	WM53
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA, 5% sucrose
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	945
Gene Symbol	CD33
Uniprot ID	CD33_HUMAN
Immunogen	Human AML cells
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081