

Anti-CD11b antibody [VIM12] {PE} (STJ16101398)
STJ16101398**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD11b is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	VIM12
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	984
Gene Symbol	CDK11B
Uniprot ID	CD11B_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	