

Anti-c-MET (Hepatocyte Growth Factor Receptor) antibody [3D4] (STJ16101402)

STJ16101402

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-c-MET (Hepatocyte Growth Factor Receptor) is suitable for use in Immunohistochemistry, Immunocytochemistry and Immunofluorescence research applications.
Applications	IHC-P/ICC/IF
Host/Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3D4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P ICC IF 1:250
Formulation	PBS with 1% BSA and 0.1% sodium azide
Isotype	IG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide derived from cytoplasmic domain of human c-Met protein
Immunogen Region	
Specificity	
Immunogen Sequence	

