

**Anti-CD120b antibody [MR2-1] (STJ16101405)**

STJ16101405

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-TNF-RII p75/p80 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	IHC-F/FC/ELISA/IF/IP/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Cynomolgus Monkey/Rhesus

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	MR2-1
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	0.2 um filtered protein G purified antibody
<b>Dilution</b>	
<b>Range</b>	
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen</b>
<b>Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>