

Anti-TCC antibody [aE11] (STJ16101406)

STJ16101406

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TCC is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunofluorescence and Immunohistochemistry research applications.
Applications	IHC-F/FC/ELISA/IF/IHC-P
Host/Source	Mouse
Reactivity	Human/Horse/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	aE11
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence