

Anti-CDK4 antibody [EP180] (STJ16101441)

STJ16101441

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-CDK4 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EP180
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:100-1:200
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	1019
Gene Symbol	CDK4
Uniprot ID	CDK4_HUMAN
Immunogen	a synthetic peptide corresponding to residues of human CDK4 protein
Immunogen Region	
Specificity	
Immunogen Sequence	