

Anti-Phospho-Akt antibody [LP18] (STJ16101675)

STJ16101675

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Phospho-Akt is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LP18
Concentration	>=300 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A chemically synthesized phosphorylated peptide corresponding to residues around Ser473 of the human Akt molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	

