

Anti-Bcl-6 Oncoprotein antibody [LN22] (STJ16101682)

STJ16101682

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Bcl-6 Oncoprotein is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LN22
Concentration	>=220 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:60
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant protein containing 1-350 amino acids of the N-terminus of the Bcl-6 human oncoprotein molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081