

Anti-Beta-Sarcoglycan antibody [betaSarc1/5B1] (STJ16101684)

STJ16101684

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Beta-Sarcoglycan is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	betaSarc1/5B1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6443
Gene Symbol	SGCB
Uniprot ID	SGCB_HUMAN
Immunogen	Fusion protein RBSG-NT of the human beta-sarcoglycan sequence.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081