

Anti-CA19-9 antibody [C241:5:1:4] (STJ16101686)

STJ16101686

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Sialyl Lewis a is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-P/IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	C241:5:1:4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:200 IHC-F
Formulation	Purified in PBS, 1% BSA and 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human colorectal adenocarcinoma cell line COLO 205.
Immunogen Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081